

ASAHI CHIKAI V

ASAHI Peripheral Vascular Guide Wire

ASAHI 0.014" Guide Wire with 3S effect

Superior torque performance

Superb shape retention

Smooth delivery




ASAHI CHIKAI V

A 0.014" Guide Wire with ACT ONE technology from ASAHI INTECC




ACT ONE
Asahi Cable Tube



Ordering Information

■ Structure



Employs the ASAHI brand rope coil, which provides torqueability, flexibility, resiliency.

Product	Catalog No.	Total Length	Outer Diameter	SLIP COAT® coating Length	Tip Shape
ASAHI CHIKAI V	WCKV165-14	165cm	0.36mm (0.014-inch)	155cm	Straight
	WCKV180-14	180cm		170cm	

Coated with SLIP-COAT® coating. SLIP-COAT® coating is protected by U.S. Pat. 5331027. SLIP-COAT® is a trademark of Angiotech BioCoatings Corp. in the U.S. and /or other countries.



Manufactured by:

Distributed by:



ASAHI INTECC CO.,LTD.
<http://www.asahi-intecc.com>