

Hemostasis Valve Adapter



Gateway HV™ Valve Adapter

The Gateway HV™ attachable hemostasis valve adapter is the answer to minimizing blood loss during the introduction of catheters or devices into the vasculature. For those times in procedures when a device needs to “park” for a period of time, The Gateway HV™ adapter is your go-to device.

The Gateway HV™ consists of a luer fit hub, a hemostasis valve, and side arm tubing with a stopcock.

Features:

- Ability to pass up to an 8Fr device
- Ability to seal on mandrel wires
- Available separate sterile
- Pressure rated up to 30 PSI

Ordering Information:

GWT-000-30 Hemostasis Valve Adapter 10/box

