


Structure & ordering information



# ASAHI SUOH 03

PTCA GUIDE WIRE

Product	Catalog No.	Diameter	Coating	Usable length	Coil length	Radiopaque length	Tip shape	Label color	Clip color
ASAHI SUOH 03	AHW14R013P	0.36mm (0.014inch)	Hydrophilic coating (SLIP-COAT®) 52cm	190cm	19cm	3cm	Pre-shape	Red	Orange and Red
	AHW14R313P	0.36mm (0.014inch)	Hydrophilic coating (SLIP-COAT®) 52cm	300cm	19cm	3cm	Pre-shape		



CE 0344

Coated with SLIP-COAT® coating.

© 2017 ASAHI INTECC CO., LTD. "ASAHI", "SUOH" and "SION" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Manufactured by:

Distributed by:



ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan  
 Contact phone number: +81-561-86-9101  
<http://www.asahi-intecc.com>



ASAHI SUOH 03/ver.1/AMC-P17016